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### Confezione refill, 1 siringa da 4,5 g:

DEI® experience Dentina:

A1 (REF. EXDA1), A2 (REF. EXDA2), A3 (REF. EXDA3), A3.5 (REF. EXDA3.5),

A4 (REF. EXDA4), B1 (REF. EXDB1), B2 (REF. EXDB2), B3 (REF. EXDB3),

B4 (REF. EXDB4), C1 (REF. EXDC1), C2 (REF. EXDC2), C3 (REF. EXDC3),

C4 (REF. EXDC4), D2 (REF. EXDD2), D3 (REF. EXDD3), D4 (REF. EXDD4).

DEI® experience Cervicale:

A2 (REF. EXCA2), A3 (REF. EXCA3), A4 (REF. EXCA4),

B4 (REF. EXCB4), C4 (REF. EXCC4), D4 (REF. EXCD4).

DEI® experience Smalto:

LIGHT (REF. EXSL), MEDIUM (REF. EXSM), DARK (REF. EXSD), IMPACT (REF. EXI 1)

DEI® experience Massa Speciale:

INCISALE BIANCO (REF. EXM50), INCISALE OPACO 90% (REF. EXM90),

INCISALE BLU (REF. EXMB), TRANSLUCIDO (REF. EXMT),

WINDOW (REF. EXMW)

DEI® experience Gengiva Flow:

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LIGHT (REF. EXG L), MEDIUM (REF. EXG M), DARK (REF. EXG D),

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Resin based on urethane dimethacrylate, barium glass filler; silica, polymer powder. Bis-GMA free.

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# DEI® experience

THE NEW COMPOSITE SYSTEM FOR  
UNIVERSAL RESTORATIONS AND  
PREPARATIONS OF BRIDGES, CROWNS,  
INLAYS, FACETS, AND CHARACTERIZATION  
OF MOBILE PROSTHESES

Bis-GMA FREE

CE  
0546

Product for dental use  
Keep away from children  
Date of the specifications: 06/2019

## CHARACTERISTICS

The planning of DEI® experience issues from the consideration that the natural tooth owns a complex chromatic characterization, due either to the shape and surface texture, or to the superimposition of different structures carrying various translucencies and opacities that must be “copied” for the purpose of restoring a correct visual perception of natural and aesthetic conditions.

The simple technique of stratification used with DEI® experience permits the achievement of ever got before aesthetic results. Peculiarity of the product is the presence of three masses having different opacities, that enable the achievement of natural, highly aesthetic restorations, using an easy and quick technique.

The presence of the mass “enamel” with its special grade of transparency, available in 4 tints, gives the tooth the natural translucence brought about by passing of light, that is reflected on the inner mass of “dentin”, available in the 16 Vita (!) colours.

The supply of 6 cervical tints more opaque than the enamels and dentins, to be placed in the neck area and used in presence of reduced depths for easier making of the colour, increases the realistic chromatic effect of the natural tooth.

## COMPOSITION

Resin based on urethane dimethacrylate, barium glass filler; silica, polymer powder.  
Bis-GMA free.

## Technical data

Resistance to compression:	350 MPa
Transverse resistance (flexion):	130 MPa
Transverse resistance modulus:	9.000 MPa
Water absorption:	20 µg/mm <sup>3</sup>
Total filler by weight:	80%

Colour stability:	optimum
Polishing:	high
Radiopaqueness:	200% AL
Cervical masses translucency:	20-35%
Dentin masses translucency:	34-46%
Enamel and incisal masses translucency:	53-62%
Particle size:	0.02 - 1 micron

## WAY OF USE (STUDIO)

Proceed as usual preparing the cavity, protecting the pulp, etching and applying the adhesive system. Use of the adhesive system New Extra Bond One (not included in the kit) is recommended.

Apply the composite material stratifying the various masses. Every layer must be cured for 40 seconds for depths not thicker than about 1.5 mm or cured for 5 sec. with DEI® experience LED Studio. Photo-polymerize with curing lamps with light performance of at least 500 mW/cm<sup>2</sup>. If it is impossible to put the light source close enough to the material, it is necessary to lengthen the exposure time. The observance of photo-polymerization times and depths are fundamental for achieving complete conversion of the materials, and ensuring resistance and duration in the time. The composites harden also very few seconds, but are not polymerized. In this phase saving time could cause failure of the restoration.

Maximum Depth of polymerization:	
for enamel mass:	3.5 mm
for cervical mass:	1.7 mm
for dentin mass:	2.5 mm
Work time at full light:	60 sec.

## ADHESION ON CURED COMPOSITE

To get perfect adhesion between already cured and finished composite and composite not yet cured (for example in case of repairs): roughen a little the already cured composite, clean it with ethylic alcohol and clean brush; apply on the so prepared surface a thin layer of DEI® experience Seal Coat Fast; don't cure and add new composite. Light cure for 60 seconds with halogen or led lamp for dentist' surgeries (or for 8 minutes with DEI® experience LED in laboratory).

## FINISHING TECHNIQUE (STUDIO)

Before the photo-polymerization of last layer of composite, apply DEI® Surface Hardener and cure for 40 seconds, wash and dry. The hardening of the composite shall be complete. Using this product avoids the necessity of eliminating the superficial layer of inhibition with burs or other products, with shorter work times.

Proceed with the modellation and with the creation of surface texture. Clean the surface

with ethilic alcohol. Apply DEI® experience Seal Coat Fast spreading it with rotational movement using microbrush or small brush. After 2 minutes remove the surplus using a new microbrush and light cure with lamp DEI® experience LED Studio for 3 minutes in dentist' surgeries (or for 8 minutes with DEI® experience LED in laboratory).

Eliminate from the surface the photo-polymerized DEI® experience Seal Coat Fast using DEI® experience Polish Paste and with goat hair brush in dentist' surgeries (or cleaner 100 µm aluminium dioxide and pumice (1:1) in laboratory). Polish with cotton bur and DEI® experience Polish Paste or Polish Felt.

## WAY OF USE (LAB)

DEI® experience is the suitable composite for the application with pressing technique, known as M.C.M. (Monolithic Composite Method). We suggest to dental technicians to join us in free Whork Shops held by the Odt. Paolo Pagliari (creator of the technique) at our headquarters.

## WARNINGS

In restorations with composite don't use products containing eugenol. Keep refrigerated at constant temperature of 4°C, or in environment with fresh and constant temperature lower than 25°C. Take the pack out from the refrigerator some hours before use. Don't expose to sources of light or heat. After use seal again the syringe with its suitable cap. To release the pressure of the material, turn the piston halfway anticlockwise. Avoid contact with the skin, as in some cases a hypersensitivity may occur: in these cases don't use. For the utilization of the accessory products, read always the instructions for use of the specific product.

Attention: the liquid for opaque contains substances that must be handled with care, and is inflammable. Don't smoke during use and don't keep close to sources of heat; avoid electrostatic charges. Keep in airy environment and with cap well sealed. In case of contact with skin, wash with plenty of water and soap; in case of contact with eyes, rinse with plenty of water and consult a specialist. The peculiar characteristics of DEI® experience allow modifications of the viscosity of the product either heating the masses for making it more fluid, or cooling down them for making it more compact.

## CONSERVATION AND VALIDITY

If correctly kept at temperature lower than 25°C, the expiration date of the product is 3 years from the date of production.

(!) trade mark registered by Vita Zahnfabrik

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