DEI® experience GENGIVA FLOW

COMPOSITO FOTOPOLIMERIZZABILE

CE
0425

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LIGHT-CURING COMPOSITE

CE
0425

Product for dental use only. Keep away from children.
Date of specifications: 04/2021

Available packages

Two 3 g syringe in colours:
LIGHT REF EXG L, MEDIUM REF EXG M,
DARK REF EXG D, SUPER DARK REF EXG SD

Areas of application

- Gingiva tissue restorations and reproductions
- Implant supra structures
- Characterization of full & partial dentures

Composition

Filler, inorganic: 57% by weight, 51% by volume
Glass filler (mean particle size): 0.7 µm
Pyrogenic silicic acid (mean particle size): 0.04 µm
Monomers (43% by weight): Urethane dimethacrylate,
Tetramethylene dimethacrylate, Bis-GMA
Additional substances (1% by weight): initiators, stabilizers, pigments

Technical data

Working times: 1-3 minutes (depending on the light conditions)
Depth of polymerization: 2 mm (with exposure for 40 sec.)
Flexural strength: 140 MPa
E-Module: 6500 MPa
Compressive strength: 370 MPa

Instructions for use [Laboratory]

The frame should be provided with retentions and treated with Dream Opaque before the application of the composite.
For the adhesion process to the Dream Opaque, check the instructions of Dream Opaque. If the opaque surface has been contaminated apply Seal Coat Fast on it, no light cure it and then the composite Gengiva Paste with a spatula or whatever you like and adapt it carefully. Light cure for 2 minutes with suitable lamp.

Instructions for use [Practice]

In case of intra-oral applications light cure the all area for 2 minutes. If the thickness of the restoration is more than 2 mm, it's important to make incremental stratifications avoiding to make higher layers that could be not correctly light cured into the deep.

Polishing

Consider the DEI® Italia protocol including the following products: Surface Hardener, Seal Coat Fast, Goat Brush, Polish Paste and Polish Felt.

Repair

In case of necessity to repair or adjust the form of light cured composite make it as follows:
1) Roughness the surface with a bur (aspire dust during the work).
2) Clean with ethylic alcohol.
3) Apply Seal Coat Fast without light curing.
4) Apply the composite and light cure.
5) Polishing (see above Polishing).

Warnings

- Avoid precontacts.
- Check that composite fits over a suitable frame.
- Apply the individual materials in the correct thickness.
- Avoid to include air bubbles during the handling of composite.
- If colour change check the efficiency and shelf life of the bulb, the respect of light curing time, thickness of the layers and the application of the correct polishing protocol.
- Close immediately after use the composite syringe.
- Do not use after the expiry date printed on the label of each syringe.
- Do not use products containing eugenol because this component will inhibits the curing.

Contra-indications

Contraindications and side effect are unknown at the moment.
If a patient has known hypersensitivities towards a component of this product, we recommend not to use it.

Storage

Store this material between 10-25°C.
If the composite is kept in a refrigerator, it must be removed well before use to resume its original contouring properties.

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