**CEMENTO COMPOSITO DUALE
ULTRA OPACO**

ISTRUZIONI D'USO

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Dimethacrylates, dental glass, silicon dioxide, photo initiators, catalysts.

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Manufacturer:
First Scientific Dental GmbH
Robert-Bosch-Strasse 26,
25335 Elmshorn, Germany

DEI®italia

Via Torino, 765 - 21020 Mercallo (VA) - Italy
Tel. +39 0331 969270 - www.deitalia.it

DUAL CURE ULTRA OPAQUE COMPOSITE CEMENT

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

DEI® Ghost Cem is an ultra opaque highly masking dual cure easy flowing microhybrid composite for permanent cementations. DEI® Ghost Cem is a **radiopaque** and relatively high filled composite of very high strength. DEI® Ghost Cem is based on methacrylate resin and inorganic filler particles of **0.05-1 µm**. The total filler load is 61% and the total filler volume 41%. DEI® Ghost Cem meets the requirements of **DIN EN ISO 4049, type 2, class 3**.

INDICATIONS / INTENDED USE

Permanent cementation of:

- Crowns, bridges and mesostructures on titanium and zirconia abutments
- Crowns and bridges on highly discolored tooth structures

PERFORMANCE FEATURES

The performance features of the product meet the requirements of the intended use.

CONTRAINDICATIONS

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

PATIENT TARGET GROUP

Persons who are treated during a dental procedure.

INTENDED USERS

This medical device should only be used by a professionally trained dental practitioner.

INCOMPATIBILITY WITH OTHER MATERIALS

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

APPLICATION

PREPARING THE MINIMIX-SYRINGE

Remove the cap of the MINIMIX-syringe and throw it away (**do not use it again!**). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The material is now ready for application. **The working time (23°C (74°F)) in the self cure mode is 2:00 minutes from start of mixing.**

NOTE

Store used syringe with fixed used mixing cannula. Discard the first 2-3 mm (about the size of a peppercorn) of the extruded material. This has to be done for each new mix.

CEMENTATION

1. PREPARING

Prepare the luting side areas of the restoration (crown, bridge or mesostructure) with a suitable bonding agent:

- Restorations made from:
- Zirconia with e.g. DEI® UniAdhesive 1
 - Silicate ceramics with e.g. DEI® Ceramic Prime
 - Metal (precious or non-precious metal) with e.g. DEI® UniAdhesive 1

2. CEMENTING OF CROWNS, BRIDGES AND MESOSTRUCTURES

a) Abutment
Prepare the abutments according to manufacturer instructions.

b) Tooth structure
Apply a suitable dual cure bonding agent to the prepared tooth surface. For cementing (preparation of the restoration see 1.) apply a uniform coating (0.5 mm) of DEI® Ghost Cem on the luting side of the restoration. Seat the restoration under light pressure and retain it in place until the cement has completely cured. Remove excess material. Light cure all marginal areas of the restoration from each direction for **20 seconds** with a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm². The dual cure cement system will auto cure within **4 minutes**.

3.1 DEI® ITALIA PROTOCOL FOR CEMENTATION

The substructures in metal, carbon or other materials must be treated as per the pre-opaque protocol, given that, with the use of this cement no pre-opacity is required. We recommend considering an internal anatomical reduction of the elements (digital design) of 0,1 mm.

3.2 DEI® ITALIA PROTOCOL FOR ELEMENTS IN ZIRCONIO, METAL ALLOYS, PMMA AND COMPOSITES

- To cement ceramic or silicate elements, replace DEI® UniAdhesive 1 with DEI® Ceramic Prime in the following protocol:
1. Sandblast at 2 atm with Biox 80/ 60 micron.
 2. Clean with ethyl alcohol and a clean brush.
 3. Once dry, apply DEI® UniAdhesive 1.
 4. Once dry, apply a very thin layer of DEI® New Extra Bond Two mixed 1:1 with DEI® New Extra Catalyst.
 5. Pass a blow of clean air (free oil and dry compressor).
 6. Light cure each element for 20 seconds with studio light and 4 minutes in a laboratory light unit for the whole structure.
 7. Cement with DEI® Ghost Cem, positioning the element on the structure, give 2/3 seconds to the excess cement then remove the excess with a probe.
 8. Light cure for 6 minutes with a laboratory light unit or 10 seconds per side (vestibolar, occlusal, internal) of the element.

STORAGE

Do not store above 20°C (68°F). Protect from direct sunlight. Store unopened material in the refrigerator. Opened cartridges have to be used up within 3 months. Do not use after expiry date.

Additional Notes / Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Color stability meets the requirements of DIN EN ISO 4049.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

COMPOSITION

Dimethacrylates, dental glass, silicon dioxide, photo initiators, catalysts

DISPOSAL

Disposal of the product according to local authority regulations.

REPORTING OBLIGATION

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

NOTE

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices. (EUDAMED - <https://ec.europa.eu/tools/eudamed>)

WARRANTY

First Scientific Dental GmbH warrants this product will be free from defects in material and manufacture. First Scientific Dental GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and First Scientific Dental GmbH's sole obligation shall be repair or replacement of the First Scientific Dental GmbH product.

LIMITATION OF LIABILITY

Except where prohibited by law, First Scientific Dental GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Manufacturer:
First Scientific Dental GmbH
Robert-Bosch-Strasse 26,
25335 Elmshorn, Germany

 **DEI® italia**

Via Torino, 765 - 21020 Mercurio (VA) - Italy
Phone +39 0331 969270 - www.deitalia.it