**DEI®
experience
seal coat
FAST**

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LIGHT CURING
SEALING COATING

Available packages:

DEI® experience Seal Coat Fast contains: one 30 ml bottle (REF. DRY 3)

Product for dental use
Keep away from children
Date of the specifications: 02/2019

Characteristics

DEI® experience Seal Coat Fast is a light curing sealing coating for composites and resins.

DEI® experience Seal Coat Fast enhances aesthetics and resistance of the treated surface and markedly reduces the growth of bacterial plaque. Tests run with atomic force microscope show that DEI® experience Seal Coat Fast is able to seal micro-porosities

and improve surface smoothness.

Possible applications:

- Restorations in composite.
- Prosthetic crowns in composite.
- Resins for pressings.
- Fastening of calcineable retentions (spread a thin layer on the crown, wait 3-4 minutes, and apply from on high the retentions; light cure for 25 seconds).

For the curing times in other types of work, follow our protocols.

DEI® experience Seal Coat Fast is also usable for repairs on composite, reactivating the bonds necessary for connecting polymerized composite and composite to use for repair.

It is also utilisable to improve the adhesion on methacrylate resins, for instance for characterizations in composite on splint teeth.

Composition

Acrylate.

Way of use in oral cavity

Before light curing the last layer of composite, apply DEI® Surface Hardener (not included) for eliminating the inhibiting surface. After the light curing, take out the gel of Surface Hardener. With a disposable brush spread a thin layer of DEI® experience Seal Coat Fast on the surface, blow well with clean and oil/water free air for better dispersion of the product. Let work 10 seconds, and light cure for 90 seconds. Remove the surface layer of DEI® experience Seal Coat Fast using DEI®

experience Polish Paste and small go-at hair brush, and polish with cotton brush (both available under request). DEI® experience Seal Coat Fast shall remain only as a seal of micro-porosity.

Extra-oral use (Lab.)

DEI® Italia Polishing protocol

The DEI® experience composite is easily polished with the help of DEI® Polish Paste. We can seal the surface making it stable and reducing plaque attachment.

In our research laboratories we have developed a protocol that makes the surface even brighter.

1. Proceed with all the controls of the MCMV handwork in the articulator.
2. Gently polish with DEI® Finishing Roll REGULAR.
3. Clean well with pure alcohol.
4. When the surface is dry, moisten the entire surface with DEI® experience Seal Coat Fast.
5. Position the element in the dark for 4 minutes.
6. Remove the excess of DEI® experience Seal Coat Fast with a dry brush or a microbrush.
7. Wait 2 minutes and light-cure for 8 minutes.
8. Polish the surface of the element with DEI® Finishing Roll FINE.
9. Apply DEI® Polish Paste with white goat toothbrush at low revolutions to avoid burning the composite.
10. Pass the cotton brush up to high revolutions but in continuous cycles to avoid burning the composite.

11. Wash with warm water and soap. A video tutorial by Paolo Pagliari is available on our website www.deitalia.it.

Use as promoter of adhesion

This protocol is suitable for composite, acrylic resin, PMMA and PEEK.

Run a bur on the surface to reactivate, clean and remove the grease (out of oral cavity ethyl alcohol can be used), dry well, run DEI® experience Seal Coat Fast. Don't light cure. Apply the composite and light cure.

Warnings

Flammable! Irritant, avoid contacts with skin, burns, eyes. In case of contact rinse with lots of water; for contact with eyes, after rinse contact immediately a specialist.

Seal well the bottle after use. Don't use in case of individual hypersensitivity to the components.

DEI® experience Seal Coat Fast contains components with different specific weights, so **the bottle must be shaken well before use**.

Storage

Keep at temperature lower than 25°C. After use seal well the bottle.

Via Torino 766 - 21020 Mercallo (VA) - Italy
Phone +39 0331 963270 - Fax +39 0331 963271
www.deitalia.it

DEI® Lab / DEI® experience

SEAL COAT FAST

RIVESTIMENTO SELLADOR
FOTOPOLIMERIZABLE

Formas de suministro:

DEI® experience Seal Coat Fast contiene: 1 frasco de 30 ml (REF. DRY 3)

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Characteristics

DEI® experience Seal Coat Fast is a re-stimulating sealator fotopolimerizable para composites and resins.

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