DEI® experience

SEAL COAT FAST LED

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Multifunctional methacrylates, ethylacetate.

DEI® experience Seal Coat Fast Led non contiene methylmethacrylate.

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DEI® experience
**SEAL COAT
FAST LED**
LIGHT CURING
SEALING COATING

Via Torino 765 - 21020 Merccallo (VA) Italy
Tel. +39 0331 969270
Fax +39 0331 969271
www.deitalia.it

DANGER

DEI® experience
**SEAL COAT
FAST LED**
RIVESTIMENTO SELLADOR
FOTOPOLIMERIZABLE

Via Torino 765 - 21020 Merccallo (VA) Italy
Tel. +39 0331 969270
Fax +39 0331 969271
www.deitalia.it

DANGER

Product for dental use
Keep away from children
Date of the specifications: 07/2019

Available packages:
DEI® experience
Seal Coat Fast Led contains:
one 5 ml bottle **[REF] DRY 2**

Characteristics
DEI® experience Seal Coat Fast Led is a light curing sealing coating for composites and resins.

DEI® experience Seal Coat Fast Led enhances aesthetics and resistance of the treated surface and markedly reduces the growth of bacterial plaque. Tests run with atomic force microscope show that DEI® experience Seal Coat Fast Led is able to seal micro-porosities and improve surface smoothness. Possible applications:

- Restorations in composite.
- Prosthetic crowns in composite.
- Resins for prostheses.
- Fastening of calcinable retentions (spread a thin layer on the crown, wait 3-4 minutes, and apply from on high the retentions; light cure for 25 seconds).

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Composition
Multifunctional methacrylates, ethylacetate.
DEI® experience Seal Coat Fast Led does not contain methylmethacrylate.

Use as promoter of adhesion
DEI® experience Seal Coat Fast Led is also usable for repairs on composite, reactivating the bonds necessary for connecting polymerized composite and composite to use for repair.

It is also utilizable to improve the adhesion on methacrylate resins, for instance for characterizations in composite on splint teeth.
PROTOCOL

Run a bur on the surface to reactivate, clean and remove the grease (out of oral cavity ethylic alcohol can be used), dry well, run DEI® experience Seal Coat Fast Led. Don't ligh cure. Apply the composite and light cure.

Way of use in oral cavity
Before light curing the last layer of composite, apply DEI® Surface Hardener (not included) for eliminating the inhibiting surface layer. After the light

Composicion
Multifunctional methacrylates, ethylacetate.
DEI® experience Seal Coat Fast Led no contiene methylmethacrylate.

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curing, take out the gel of Surface Hardener. With a disposable brush spread a thin layer of DEI® experience Seal Coat Fast on the surface, blow with clean and oil/water free air for better dispersion of the product. Let work 10 seconds, and light cure for 90 seconds. Remove the surface layer of DEI® experience Polish Paste and small go at hair brush, and polish with cotton brush (both available under request). DEI® experience Seal Coat Fast shall remain only as a seal of micro-porosity.

Warnings
Flammable! Irritant, avoid contacts with skin, gums, eyes. In case of contact rinse with plenty of water; for contact with eyes, after rinse contact immediately a specialist.

Seal well the bottle after use. Don't use in case of individual hypersensitivity to the components.

Storage
Keep at temperature lower than 25°C. After use seal well the bottle.

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