# DEI® CLEVER REPLY NEW

## COMPOSITO PER RESTAURI ESTETICI

CE  
0425

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DEI® Clever Reply New:

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# DEI® CLEVER REPLY NEW COMPOSITE FOR AESTHETIC RESTORATIONS

Date of the specifications: 04/2021

## AVAILABLE PACKAGES

DEI® Clever Reply New:

One 4,5 g syringe in the following colours:  
A1 (REF. COM 70), A2 (REF. COM 71), A3 (REF. COM 72),  
A3,5 (REF. COM 73), A4 (REF. COM 85), B2 (REF. COM 75),  
B3 (REF. COM 76), C2 (REF. COM 77), C3 (REF. COM 78),  
D3 (REF. COM 79).

## CHARACTERISTICS AND CLASSIFICATION

DEI® Clever Reply New is the new DEI® italia hybrid composite developed for attaining the even higher aesthetic performances wished by dentists who consider the aesthetic conservation to be a reconstructive art.

Therefore, DEI® Clever Reply New offers more possibilities of choice of colours with 10 tints in the introductory kit. DEI® Clever Reply New is recommended for I, II, III, IV V classes of restorations for reconstructions of facets and stumps.

DEI® Clever Reply New does not stick to instruments, is coating, therefore guarantees best working. We recommend to use through the various work steps of restoration DEI® New Extra Bond One and the finishers (not included) DEI® Surface Hardener and DEI® Seal Coat Fast.

## COMPOSITION

Glass powder, diurethane dimethacrylate, silicon dioxide, Bis-GMA, tetramethylene dimethacrylate.  
Total filler: 75% by weight (53% by volume) anorganic filler (0.005 - 3.0 µm)

## TECHNICAL DATA

Vickers hardness: ≥ 539 MPa  
Flexural strength: ≥ 120 MPa  
Modulus of elasticity: ≥ 9000 MPa  
Compressive strength: ≥ 400 MPa  
Water absorption: ≤ 40 µg/mm<sup>3</sup>  
Water solubility: ≤ 7.7 µg/mm<sup>3</sup>  
Depth of polymerization: ≥ 1.5 mm  
Filler content by weight: 75%  
Filler content by volume: 53%  
Radiopacity: 100% AI  
Particle size distribution: 0.005 to 3.0 µm

## WAY OF UTILIZATION

### Preparation of the cavity

Isolate with dam the operating area. The preparation must not be contaminated by saliva or other items. Prepare the cavity seeking to retain as much enamel as possible. To improve the closing of the rims and the final safety of the restoration, chamfer at 45-60 degrees the edges. Preparing rear teeth, avoid to leave enamel without support of dentin.

### Technique of etching

Apply the etching acid only on the enamel, taking care to cover the chamfer and not wet the dentin. Let the acid work for 30 seconds. Rinse for 20 seconds. Dry with air. The etched enamel shall have a chalky aspect; otherwise, or in case of contamination by oil, water or saliva, repeat the process. The etching acid is irritant for skin, mucosas and eyes. Consequently, avoid contacts. Eventually, rinse with plenty of water and, if necessary, call for a specialist.

### Technique for use of DEI® New Extra Bond One

Apply DEI® New Extra Bond One or Two following the instructions.

### Application of composite

Apply layers of composite not thicker than 1.5 mm. Light cure for 40 seconds (or 10 seconds with DEI® experience LED Studio). Repeat the procedure up to full restoration of the tooth.

## Finishing technique

Before light curing the last layer of composite, apply DEI® Surface Hardener and light cure for 40 seconds, wash and dry. The hardening of the composite shall be complete. The use of that product avoids the need of eliminating the surface inhibition layer with burs or other items, consequently with shorter working times.  
After modeling, run a micro-brush with ethylic alcohol on the surface. With disposable brush spread a thin layer of DEI® Seal Coat Fast Led on the surface, blow with clean and oil/water free air to improve the distribution of the product; run a clean brush or micro-brush for elimination of surpluses, let work 10 seconds and light cure for 90 seconds. Take out the surface layer of DEI® Seal Coat Fast Led using DEI® Polish Paste and goat hair brush, and polish with cotton small brush. DEI® Seal Coat Fast Led remains sealing the micro-porosity.  
At the end of polishing an extremely hard, resistant and shiny surface shall be obtained. In addition, the present technique allows considerable diminution of growth of tartar, bacterial plaque and of colorations caused by beverages such as coffee.

## SPECIAL NOTES

- The working time under a surgical lamp is approximately 2 minutes.
- In case of time-consuming restorations, the surgical lamp should be either temporarily moved away from the working area or the material should be covered by an opaque foil in order to prevent the composite from curing too early.
- Use a light-curing unit with an emission spectrum of 350 - 500 nm (Light intensity for curing: ≥ 650 mW/cm<sup>2</sup>) for the polymerization of this material. As the required physical properties can only be achieved when the lamp is not obsolete, works correctly and its luminous intensity must be checked regularly as described by the manufacturer.

## HAZARD AND PRECAUTIONARY STATEMENTS

### Contains tetramethylene dimethacrylate

**Warning:** May cause an allergic skin reaction. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

## CONTRAINDICATIONS / INTERACTIONS

If a patient has known allergies against or hypersensitivities to a component of this product, it may not be used or only under strict medical supervision by the doctor/dentist. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

## SIDE-EFFECTS

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product - please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentin is exposed, the pulp must be protected adequately (e. g. calcium hydroxide preparation).

## INTERACTIONS WITH OTHER SUBSTANCES

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zinc-oxide eugenol cements) containing such substances.

## STORAGE

Store at 10-25°C (50-77°F). Close the screw syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures from becoming blocked.

Do not use after expiry date (refer to label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

\* Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

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PROBLEM	CAUSE	REMEDY
Composite does not cure properly	Light output of the light-curing lamp is inadequate Emitted wavelength range of the lightcuring lamp is inadequate	Check the light output and change the light source if required Consult the manufacturer of the lightcuring lamp. Recommended wavelength range: 350-500 nm
Composite in the syringe is sticky and soft, colorless liquid separates in the syringe	Material has been stored for a longer period at 25°C (77°F) Material has been kept in a syringe warmer for too much time	Adhere to storage temperature. Store at 10-25°C (50-77°F) Never keep a syringe in a syringe warmer for more than one hour per application
Composite appears too hard and firm in the syringe	Material stored at temperatures, 10°C (50°F) for a longer period of time Syringe not properly sealed, composite partially cured	Allow the composite to heat to room temperature before use; use a syringe warmer if necessary Always seal the syringe properly with the cap after taking out composite
Inlay/onlay is not properly retained when fitted	Restoration is too opaque to be cemented using only light-curing composite	Use dual-curing luting composite
Composite does not cure completely (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. thickness of 2.0 mm per layer
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 sec.

Via Torino 765 - 21020 Mercurio (VA) - Italy  
Phone +39 0331 969270 - Fax +39 0331 969271  
www.deitalia.it

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